



LNP 2.0 platform for RNA Technology

# PRODUCT DATA SHEET

Document No.: PDS-NV1012    Revision: 1.0    Effective Date: JUL/2026

## 1 Product Identification

Company Logo &  
Contact info



info@neovac.co.uk

Product Name:            **LipTx™ NV1-012**

Catalog Number:        **NV1012**

## 2 Product Description

### 2.1 Short summary

**LipTx™ NV1-012** is a proprietary ionizable lipid developed by NeoVac for nucleic acid delivery applications. The lipid is designed to facilitate efficient nucleic acid encapsulation during LNP manufacturing while supporting intracellular delivery through pH-dependent ionization.

### 2.2 Intended use

For Research Use Only (RUO).

Use for Lipid Nanoparticle (LNP) Formulation for nucleic acid delivery.

Not intended for any human use, therapeutic administration, diagnostic procedures, food, or veterinary use.

### 2.3 Availability and lead times

| Grade                   | Availability           |
|-------------------------|------------------------|
| Research Use Only (RUO) | On Stock               |
| GMP Grade Material      | Available Upon Request |



## PRODUCT DATA SHEET

Document No.: PDS-NV1012    Revision: 1.0    Effective Date: JUL/2026

| Catalog Number | Quantity | Price [USD] |
|----------------|----------|-------------|
| NV1012-025MG   | 25 mg    | 650         |
| NV1012-100MG   | 100 mg   | 1,950       |
| NV1012-001G    | 1 g      | 9,950       |

Lead times:

- Gram Production: 6 – 8 weeks
- GMP Grade Material: 6 -12 months

### 3 Product Specifications

#### 3.1 Molecular weight

824.37 g/mol

#### 3.2 Purity

Purity ≥95 %.

#### 3.3 Appearance

Pale yellow colour viscous liquid.

#### 3.4 Solubility profile

Max. 50 mg/ml in: Ethanol (EtOH), Methanol (MeOH), Chloroform, Isopropanol (IPA) and Dimethyl sulfoxide (DMSO).

### 4 Storage and Stability

#### 4.1 Storage

Store at –80°C in tightly closed container for best shelf life.



LNP 2.0 platform for RNA Technology

# PRODUCT DATA SHEET

Document No.: PDS-NV1012   Revision: 1.0   Effective Date: JUL/2026

For storage at  $-20^{\circ}\text{C}$  more frequent assay and impurities retests are recommended.

Recommended:

- Store in original container.
- Protect from moisture and oxygen.
- Minimize freeze-thaw cycles.
- Protect from prolonged exposure to light.

## 4.2 Shelf life

Initial shelf life at  $-80^{\circ}\text{C}$  12 months. Retests of assay and impurities for shelf-life extension possible.

## 5 Troubleshooting Guide

### 5.1 Documentation

- Material Safety Data Sheet (MSDS)
- Certificate of Analysis (CoA)
- Regulatory Package (Optional) – available under the CDA.

### 5.2 Legal and Ordering information

This product is supplied for research and development use only. **NOT FOR CLINICAL USE.**

### 5.3 Disclaimer

NeoVac products are supplied For Research Use Only (RUO) and are not intended for use in humans or animals for therapeutic, diagnostic, or clinical purposes.

These products are intended solely for laboratory research by qualified personnel using appropriate laboratory practices and in accordance with the applicable Instructions for Use (IFU), Material Safety Data Sheet (MSDS), and all relevant product documentation. Users are solely responsible for determining the suitability of the product for their intended research application and for ensuring compliance with all applicable laws, regulations, and institutional requirements.





LNP 2.0 platform for RNA Technology

## PRODUCT DATA SHEET

Document No.: PDS-NV1012   Revision: 1.0   Effective Date: JUL/2026

NeoVac shall not be liable for any direct, indirect, incidental, consequential, or special damages arising from the misuse of the product, use outside its intended purpose, unauthorized modification or alteration, improper storage or handling, or failure to follow the applicable Instructions for Use (IFU) or other product documentation. Any use of the product that deviates from the manufacturer's published specifications or recommendations is undertaken solely at the user's own risk.

Except as expressly stated in the applicable product documentation or required by applicable law, NeoVac makes no representations or warranties, express or implied, including, without limitation, any implied warranties of merchantability, fitness for a particular purpose, or non-infringement.

5.4    Contact and Technical support

[info@neovac.co.uk](mailto:info@neovac.co.uk)

